

Simulator Checklist for MD 128

Central Venous Catheter Insertion Training Simulator

Note: Make sure all the following items listed below are in the carrying case/container before re-shelving the simulator

- **1 Upper Male Torso (with detachable head and lower arm pieces)**
- **1 Landmark Puncture Pad (with connected tubes--inserted into torso)**
- **1 Ultrasound Puncture Pad (with connected tubes)**
- **1 Transparent Cannulation Training Block (skin covered white base which shows veins in red and green colors —stored underneath removable beige colored block to protect it)**
- **1 big beige colored storage block to keep torso and other accessories in**
- **1 red pulsation bulb connected to tube (in plastic bag)**
- **3 plastic bottles**

(continued on back side)

- 2 blood powder bottles (1 red, 1 blue)
- 1 bottle of Sterile Saline Solution
- 1 needle syringe for fluid supply (inserted into plastic sleeve)
- Set of Two 30cm long Central Venous Catheters (kept in plastic tray in plastic bag)
- Set of 2 thin tubes (in bag with pulsation bulb)
- 1 Spoon (silver)
- Pack of large white moisture wipes
- 1 white pillow to protect simulator
- 1 large white foam pad

The Following Items Should Be Inside White Envelope

- 1 twenty-page CVC Instruction Manual
- 1 fifteen-page “Intro for CVC Training manual”
- 1 laminated “Please Read the Manual Carefully Before Use” sheet
- 1 laminated “Instruction Manual: Ultrasound Puncture Pad for M93U” sheet
- 1 three-page “Medical Definitions for MD 128” list
- 1 “Avoid Damage to Simulator” insert

